

Respireo

We are all different. The **Respireo Primo** series is a range of adjustable masks specially designed to give a personalized fit to your patients and adapt to each facial shape.

Simple and comfortable, **Respireo Primo N** is an adjustable nasal mask that will successfully assist your patients equipped with CPAP (Continuous Positive Airway Pressure) or NIV (Non Invasive Ventilation) ventilators.

Respireo Primo N is also available without the comfort adjustment dial: which is the **Respireo N**. Its silence and simplicity will help your patients to comply to their treatment every day.

In addition, thanks to their Italian fabrication, **Respireo** masks are made of high quality materials that will ensure a comfortable and long-lasting usage.



KM215000 > RESPIREO PRIMO N - S

KM215100 > RESPIREO PRIMO N - M

KM215200 > RESPIREO PRIMO N - L

KM213000 > RESPIREO N - S

KM213100 > RESPIREO N - M

KM213200 > RESPIREO N - L

KM217700 > BLUE BENDS NON VENTED X10

For further information, please contact your sales representative.



Contact :

Air Liquide Medical Systems
Parc de Haute Technologie
6, rue Georges Besse
92182 Antony cedex - France
Tél : +33 (0)1.40.96.66.00
Fax : +33(0)1.40.96.67.00

Respireo

Primo **N**

Simple and comfortable



This document is not contractual and may be amended at any time / YN091212 - 1.1 - 06/07/2011 - Raffaello Fischer / Photo credit : Le Square, Thinkstock

Respireo Primo N

1 The well-tested, robust **forehead support** ensures a personalized fit while avoiding nasal pressure for **complete patient satisfaction**.

Forehead support



2 A **dual-wall silicon** cushion preserves the security and **seal** of the mask whilst ensuring **the best comfort available** for the patient.

Dual-wall cushion



3 The attachment **hook** system makes applying the mask easier with fast removal for **minimized patient effort**.



Attachment hook system



A large and resistant **headgear** disperses pressure points for **comfort beyond the mask itself**.



Comfort dial



4 A **comfort dial** for mask adjustment to the patient allows unrivalled **precision**.



5 The mask comprises **only 2 parts** that make dismantling and assembling the mask easy, as well as simple to clean.


Dismantling and cleaning



6 Air **vent** is designed to direct a **laminar flow** downwards to **maximize quietness** during the night.



Laminar air flow vent



Comfortable with its dual-wall cushion, **Respireo Primo N** comes with an attachment system that makes it very **convenient** to use.

It is appreciated for its **silence** and will allow your patients to spend "problem free" nights with their treatment.